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 Site Distribution: SB Copy Location:	Department of Microbiology  <b>Microbiology</b>	MIC-5.5.5.3.1 Version: 1.4 Page 1 of 2 Effective: 2012/08/20  E-Authorized by: Medical Chief, Microbiology
--	---	--

### **Reporting Possible Contaminants in Blood Cultures**

Possible contaminants may include coagulase negative staphylococci, viridans streptococci, aerobic spore-forming bacilli, *Corynebacterium* species, *Propionibacterium* species, and *Micrococcus* species.

**Note:** For Paediatric blood cultures *Bacillus cereus* must be ruled out when an aerobic spore-forming bacilli is isolated.

### **Result Reporting**

#### **1. One set received:**

For coagulase negative staphylococcus and viridans streptococcus :

.....isolated after....days. Only one blood culture set received.  
Possible contaminant. Further work up will be done on specific request.  
(Text Code: **ONED**)

For other organisms :

.....isolated after....days. Only one blood culture set received.  
Possible contaminant.  
(Text Code: **ONEC**)

#### **2. Two or more strains from one set only (based on colonial morphology):**

For all organisms:

.....isolated after....days. 2 types. Possible contaminants.  
(Text Code: **TWOC**)

#### **3. Two or more sets received (One set positive out of two or three):**

For all organisms:

.....isolated after....days from one blood culture set only. Possible contaminant.  
(Text Code: **ONEO- PSBC**)

*Uncontrolled copy valid only September 9, 2015*

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 Site Distribution: SB Copy Location:	Department of Microbiology  <b>Microbiology</b>	MIC-5.5.5.3.1 Version: 1.4 Page 2 of 2 Effective: 2012/08/20  E-Authorized by: Medical Chief, Microbiology
--	---	--

### **Reporting Possible Contaminants in Blood Cultures**

#### **4. Different strains from each of two or more sets:**

##### For all organisms:

.....isolated after....days. Possible contaminant(s). This isolate is different from accession # .....  
(Text Code: **PSBC - ISODF**)

#### **5. Paediatric blood cultures**

##### For all organisms:

.....isolated after.....days. Possible contaminant(s). Further work up will be done only on specific request.  
(Test Code: **ONEP**)

*Uncontrolled copy valid only September 9, 2015*

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 Site Distribution: SB Copy Location:	Department of Microbiology  <b>Microbiology</b>	MIC-5.5.15.1.1.0 Version:13.0 Page 1 of 2 Effective: 2015/06/11  E-Authorized by: Medical Chief, Microbiology
--	---	---

### **Antibiotic Susceptibility Testing**

#### **Glossary of Codes**

AMPC: This organism is resistant to all penicillins, beta-lactamase inhibitors and all cephalosporins.

ASTND: Unable to perform antibiotic susceptibility tests due to poor growth on all AST Media.

BHSP: Beta-haemolytic streptococci are empirically susceptible to Penicillin.

BORSA: This strain shows borderline oxacillin resistance. If advice on therapy is required, please contact the microbiologist.

CREP2: This organism is resistant to all penicillins, beta-lactamase inhibitors, all cephalosporins and carbapenems.

CORAST: Interpretations are based on achievable blood levels, not on achievable ocular levels.

CPO: This organism produces carbapenemase.

CSYNR: Gentamicin resistant result indicates it will not provide synergy with a Penicillin.

ENTFLG: The majority of Enterococcus faecalis are usually susceptible to Ampicillin.

ENTFLU: The majority of Enterococcus faecalis are usually susceptible to Ampicillin and Nitrofurantoin.

Free text: For Enterococcus species, other than E. faecalis and E. faecium:

Urines: Enterococcus species are usually susceptible to Ampicillin and Nitrofurantoin.

Others (including neonatal urines): Enterococcus species are usually susceptible to Ampicillin.

ESBLP2: This organism is resistant to all penicillins, beta-lactamase inhibitors and all cephalosporins.

GNARO: Resistance to all penicillins, beta-lactamase inhibitors and cephalosporins may develop during therapy with these agents.

IRFLU: This organism is intrinsically resistant to Fluconazole.

*Uncontrolled copy valid only on September 9, 2015.*

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 Site Distribution: SB Copy Location:	Department of Microbiology <b>Microbiology</b>	MIC-5.5.15.1.1.0 Version:13.0 Page 2 of 2 Effective: 2015/06/11 E-Authorized by: Medical Chief, Microbiology
--	---	---

### **Antibiotic Susceptibility Testing**

LISTS: *Listeria* species are predictably susceptible to Ampicillin and Penicillin.

MOR2: *Moraxella* species are empirically susceptible to Penicillin.

MSEN: *Moraxella* (*Branhamella*) *catarrhalis* is predictably susceptible to cephalosporins, sulphamethoxazole/trimethoprim, erythromycin and ciprofloxacin.

NOCLSI: No clinical data available to assist in interpretation of MIC results. There is no current CLSI standard testing protocol for this organism and drug combination.

NOIVS: Susceptibility testing for this organism is unreliable. If advice on therapy is required contact the microbiologist.

PCSP: *Pasteurella canis* is empirically susceptible to Ampicillin and Penicillin

PMSP: *Pasteurella multocida* is empirically susceptible to Ampicillin and Penicillin.

SSAPU: *Staphylococcus saprophyticus* is empirically susceptible to urinary concentrations of Nitrofurantoin, sulphamethoxazole/trimethoprim and fluoroquinolones.

STANP: *Streptococcus anginosus* group are usually susceptible to Penicillin.

SUSM: Susceptibility interpretation based on meningeal infection.

SUSNM: Susceptibility interpretation based on non-meningeal infection.

*Uncontrolled copy valid only on September 9, 2015.*

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COLLECTION DATE: 17/10/15 TIME: 0620	SPECIMEN #:
ORDERED: BLOOD CULTURE SOURCE: BLOOD	
COMMENTS: Collection Site VENOUS	
BLOOD CULTURE SPECIMEN Final	MSH
*NO GROWTH DETECTED AFTER 48 HOURS INCUBATION: FURTHER REPORT TO FOLLOW ONLY IF POSITIVE*	

COLLECTION DATE: 10/10/15 TIME: 1155	SPECIMEN #:
ORDERED: URINE CULTURE SOURCE: URINE CATHETER	
URINE SPECIMEN Final CULTURE	SNY
10 to 100 X E6 CFU/L Klebsiella oxytoca 10 to 100 X E6 CFU/L Enterococcus faecalis . COMMENT: The majority of Enterococcus faecalis are usually susceptible to Ampicillin and Nitrofurantoin.	
Organism 1	KLEBSIELLA OXYTOCA
KLEBSIELLA OXYTOCA	REACTION
AMPICILLIN	R
CEFAZOLIN	I
CIPROFLOXACIN	S
NITROFURANTOIN	S
GENTAMICIN	S
TRIMETHOPRIM/SULPHAMETHOXAZOLE	S

COLLECTION DATE: 13/10/15 TIME: 1306	SPECIMEN #:
ORDERED: BLOOD CULTURE BLOOD CULTURE SPECIMEN TO SDL SOURCE: BLOOD	
BLOOD CULTURE SPECIMEN Final MSH GRAM STAIN	MSH
GRAM POSITIVE COCCI IN GROUPS SEEN IN ANAEROBIC BOTTLE INCUBATED 1.0 DAY CULTURE TO FOLLOW	
BLOOD CULTURE SPECIMEN TO SDL Final Organism 1	SNY
CULTURE	COAGULASE NEGAT. STAPHYLOCOCCI
Coagulase negative staphylococci isolated after 1 day from one blood culture set only. Possible contaminant.	

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	ORDERED: BLOOD CULTURE BLOOD CULTURE SPECIMEN TO SDL SOURCE: BLOOD	
	BLOOD CULTURE SPECIMEN Final MSH GRAM STAIN	MSH GRAM NEGATIVE BACILLI SEEN IN ANAEROBIC BOTTLE AND GRAM POSITIVE COCCI IN GROUPS SEEN IN AEROBIC BOTTLE AFTER 0.7 DAY INCUBATION CULTURE TO FOLLOW  ***** * This is a corrected result. * ***** A prior result that was reported as final has been changed. ***** * This is a corrected result. * ***** A prior result that was reported as final has been changed.  BLOOD CULTURE SPECIMEN TO SDL Final Organism 1 ESCHERCHIA COLI-ESBL PRODUCER Organism 2 COAGULASE NEGAT. STAPHYLOCOCCI  CULTURE <div style="border: 2px solid red; padding: 2px;">           Escherichia coli (ESBL producer)            This organism is resistant to all            penicillins, beta lactamase            inhibitors and all cephalosporins.         </div> isolated after < 24 hours Coagulase negative staphylococci isolated after < 24 hours Possible contaminant.

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©Imprimeur de la Reine pour l'Ontario, 2016